

The GT260 TreadCare® 'Power Walker' Rehabilitation Treadmill for home use, specifically designed for use by the elderly and patients undergoing rehabilitation. Unlike other rehab treadmills available on the market, the GT260 has a start-up speed of a just 0.1Km/h adjustable in speed increments of 0.1km/h. This design feature makes way for a very safe user platform.

The full-length side handrails are located 650mm apart (and 820mm in height from the running deck), provide extra user support and safety during use.

The GT260 ensures a safe space to walk on - with its large walking area measuring 380mm x 1070mm, this machine accommodates users weighing up to 120Kg.

To sum things up, this is one good all-round space saving treadmill ideal for home rehabilitation use.

**TREADCARE® GT260 SPECIFICATIONS:**

- Safety Stop Button and Safety Stop Tether for Emergency Power Down
- Walking area: 380mm x 1070mm
- Maximum User Weight: 120kg
- Speed range: 0.1km/h to 6km/h
- Motor Rating: 1.75HP
- Display Console/On Board Computer: Speed, Time, Distance, Calories and User Program
- Product Weight: 42Kg
- Mains Powered
- CE Certified
- Unit Dimensions (LxWxH): 1460x650x1160 (mm)
- **FREE OPTION: Delivered pre-Installed**



**BODY-SOLID UK • NORTH WEALD, ESSEX CM16 6AA**

**Tel: 0871 200 2477**

---

**[www.bodysolid.co.uk](http://www.bodysolid.co.uk)**